



Balanced #5 Forceps - Standard Tips/Straight/15cm



The Balanced #5 Forceps of the Exclusive Line collection feature an ergonomic design and ultra fine tips for optimal comfort and accuracy. The balanced action of the forceps ensures seamless efficiency during prolonged tasks, due to its rounded handle and counterweight technology providing enhanced control and precision within fine adjustments.

General Info

|                |                                               |
|----------------|-----------------------------------------------|
| Product Number | 11851-20                                      |
| Category       | Exclusive Line<br>Forceps & Tweezers > Smooth |
| EAN            | 4251909511951                                 |



Find 11851-20 on  
[finescience.com](https://finescience.com)

Technical Info

|                        |                                 |                 |        |
|------------------------|---------------------------------|-----------------|--------|
| Shipping Weight        | 65 g                            | Net Weight      | 32.5 g |
| Package Size           | 1 Pc.                           | Unit of Measure | EA     |
| Minimum Order Quantity | 1                               |                 |        |
| Alloymaterial          | Surgical Stainless Steel (Inox) |                 |        |
| Autoclavesafe          | Yes                             |                 |        |
| Length                 | 15 cm                           |                 |        |
| Style                  | Balanced #5                     |                 |        |
| Tipdimensions          | 0.1 x 0.06 mm                   |                 |        |
| Tips                   | Standard                        |                 |        |
| Tipshape               | Straight                        |                 |        |