



The Balanced #5 Forceps of the Exclusive Line collection feature an ergonomic design and ultra fine tips for optimal comfort and accuracy. The balanced action of the forceps ensures seamless efficiency during prolonged tasks, due to its rounded handle and counterweight technology providing enhanced control and precision within fine adjustments.

General Info

Product Number	11852-20
Category	Exclusive Line Forceps & Tweezers > Smooth
EAN	4251909511968



Find 11852-20 on
finescience.com

Technical Info

Shipping Weight	65 g	Net Weight	32.5 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Length	15 cm		
Style	Balanced #5		
Tipdimensions	0.05 x 0.02 mm		
Tips	Biology		
Tipshape	Straight		