



#5/30 Gold Plated Forceps - 30° Angled/Gold Plated Stainless Steel/11.5 cm



Gold-Plating gives corrosion resistance to most chemicals, salts and acids. These forceps are perfect for any application where corrosion resistance and toughness are primary requirements.

General Info

Product Number	11911-11
Category	Forceps & Tweezers > Smooth
EAN	4251909510121



Find 11911-11 on finescience.com

Technical Info

Shipping Weight	20 g	Net Weight	12 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Gold coated Surgical Stainless Steel (Inox)		
Length	11 cm		
Tipdimensions	0.1 x 0.09 mm		
Tips	Standard		
Tipshape	Angled 30°		