



#7 Gold Plated Forceps - Curved/Gold Plated Stainless Steel/12 cm



Gold-Plating gives corrosion resistance to most chemicals, salts and acids. These forceps are perfect for any application where corrosion resistance and toughness are primary requirements.

General Info

Product Number	11917-12
Category	Forceps & Tweezers > Smooth
EAN	4251909510114



Find 11917-12 on
finescience.com

Technical Info

Shipping Weight	15 g	Net Weight	15 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Gold coated Surgical Stainless Steel (Inox)		
Length	12 cm		
Tipdimensions	0.2 x 0.1 mm		
Tips	Standard		
Tipshape	Curved		