



#27 Mirror Finish Forceps - Standard Tips/Straight/13cm



These forceps have long and narrow shanks with thin, highly polished tips to prevent tissue from sticking. Ideal for embryos, nerves, or other adhesive tissues. All featured products are limited in quantity and may not be available in the future or part of the regular instrument line. While supplies last.

General Info

Product Number	11923-13
Category	Forceps & Tweezers > Smooth
EAN	4251909510138



Find 11923-13 on
finescience.com

Technical Info

Shipping Weight	15 g	Net Weight	15 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	13 cm		
Tipdimensions	0.3 x 0.15 mm		
Tips	Standard		
Tipshape	Straight		