




Mini #3E Forceps - Biology Tips/Straight/8cm/Superalloy



Mini tweezers are suitable for delicate applications which require holding a very small component firmly while granting good visibility. The short body allows a firmer grip combined with greater practicality. The superalloy material provides improved performance in severe environments.

General Info

Product Number	11961-08
Category	Forceps & Tweezers > Smooth MedTech Line
EAN	4251909512279



Find 11961-08 on
finescience.com

Technical Info

Shipping Weight	9 g	Net Weight	6 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloy material	Anti-Magnetic Anti-Acid Superalloy (CX), polished		
Antimagnetic	Yes		
Length	80 mm		
Style	#3E		
Tip dimensions	0.1 x 0.04 mm		
Tips	Standard		
Tip shape	Straight		