



Mini #2 Forceps - Standard Tips/Straight/9cm/Stainless Steel



Mini tweezers are suitable for delicate applications which require holding a small component firmly while granting good visibility. The short body allows firmer grip combined with greater practicality.

General Info

Product Number	11963-09
Category	Forceps & Tweezers > Smooth MedTech Line
EAN	4251909512293



Find 11963-09 on
finescience.com

Technical Info

Shipping Weight	11 g	Net Weight	8 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Anti-magnetic anti-acid Surgical Stainless Steel (Inox) (sa), sandblasted		
Antimagnetic	Yes		
Autoclavesafe	Yes		
Length	90 mm		
Style	#2		
Tipdimensions	0.2 x 0.15 mm		
Tips	Standard		
Tipshape	Straight		