



Mini #2A Forceps - Standard Tips/Straight/7cm/Stainless Steel



Mini tweezers are suitable for delicate applications which require holding a flat component firmly. The short body allows a firmer grip combined with greater practicality.

General Info

Product Number	11964-07
Category	Forceps & Tweezers > Smooth MedTech Line
EAN	4251909512309



Find 11964-07 on finescience.com

Technical Info

Shipping Weight	9 g	Net Weight	6 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Anti-magnetic anti-acid Surgical Stainless Steel (Inox) (sa), sandblasted		
Antimagnetic	Yes		
Autoclavesafe	Yes		
Length	70 mm		
Style	#2A		
Tipdimensions	1.23 x 0.09 mm		
Tips	Standard		
Tipshape	Straight		