



Forceps Kit - Plastic/Four tip shapes/11.5 cm



Kit includes four different full plastic polyvinylidene fluoride (PVDF) forceps. Each forcep has a length of 11.5 cm and is electrostatic discharge (ESD) safe. The four tip shapes are as follows: Straight, flat, round- Straight, fine, shar p- Curved, angled, pointed- Angled, squared, flat

General Info

Product Number	11970-90
Category	Forceps & Tweezers > Smooth
EAN	4251909510152



Find 11970-90 on finescience.com

Technical Info

Shipping Weight	40 g	Net Weight	37 g
Package Size	1 Kit	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	PVDF		
Length	11.5 mm		
Tipdimensions	Various		
Tips	Smooth		
Tipshape	Straight & Curved		