



#259 Replaceable Tips Forceps - Straight/0.5 x 0.6 mm/13 cm



Stainless steel body with replaceable plastic tips. Tips are made of polyvinylidene fluoride (PVDF) and are electrostatic discharge (ESD) safe. Body and plastic tips are corrosion and abrasion resistant. Screwdriver included. All featured products are limited in quantity and may not be available in the future or part of the regular instrument line.

General Info

|                |                                          |
|----------------|------------------------------------------|
| Product Number | 11980-13                                 |
| Category       | Forceps & Tweezers > Smooth MedTech Line |
| EAN            | 4251909510169                            |



Find 11980-13 on [finescience.com](https://finescience.com)

Technical Info

|                        |                                                        |                 |      |
|------------------------|--------------------------------------------------------|-----------------|------|
| Shipping Weight        | 30 g                                                   | Net Weight      | 30 g |
| Package Size           | 1 Pc.                                                  | Unit of Measure | EA   |
| Minimum Order Quantity | 1                                                      |                 |      |
| Alloymaterial          | Surgical Stainless Steel (Inox) body with plastic tips |                 |      |
| Length                 | 13 cm                                                  |                 |      |
| Tipdimensions          | 0.5 x 0.6 mm                                           |                 |      |
| Tips                   | Replaceable                                            |                 |      |
| Tipshape               | Straight                                               |                 |      |