



#259 Plastic Replacement Tips - 0.6 x 0.6 mm/4cm/Straight/Strong Pointed

Pack of (2) PVDF tips and (3) screws. Tips are made of polyvinylidene fluoride (PVDF) and are electrostatic discharge (ESD) safe. Compatible with Item No. 11980-13



General Info

Product Number	11988-01
Category	Forceps & Tweezers > Smooth MedTech Line
EAN	4251909510176



Find 11988-01 on finescience.com

Technical Info

Shipping Weight	5 g	Net Weight	3 g
Package Size	1 Pk	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	PVDF		
Antimagnetic	Yes		
Length	4 cm		
Tipdimensions	0.6 x 0.6 mm		
Tips	Smooth		
Tipshape	Straight / Fine		