



#5 Plastic Replacement Tips - 0.6 x 0.45 mm/4cm/Straight/Very Fine



Pack of (2) PVDF tips and (3) screws. Tips are made of polyvinylidene fluoride (PVDF) and are electrostatic discharge (ESD) safe. Compatible with Item No. 11980-13

General Info

Product Number	11988-02
Category	Forceps & Tweezers > Smooth MedTech Line
EAN	4251909510183



Find 11988-02 on
finescience.com

Technical Info

Shipping Weight	5 g	Net Weight	3 g
Package Size	1 Pk	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	PVDF		
Antimagnetic	Yes		
Length	4 cm		
Tipdimensions	0.6 x 0.45 mm		
Tips	Smooth		
Tipshape	Straight / Very fine		