



#7 Plastic Replacement Tips - 0.5 x 0.5 mm/4cm/Curved/Very Fine



Pack of (2) PVDF tips and (3) screws. Tips are made of polyvinylidene fluoride (PVDF) and are electrostatic discharge (ESD) safe. Compatible with Item No. 11980-13

General Info

Product Number	11988-03
Category	Forceps & Tweezers > Smooth
EAN	4251909510190



Find 11988-03 on
finescience.com

Technical Info

Shipping Weight	5 g	Net Weight	3 g
Package Size	1 Pk	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	PVDF		
Antimagnetic	Yes		
Length	4 cm		
Tipdimensions	0.5 x 0.5 mm		
Tips	Smooth		
Tipshape	Curved / Very fine		