



#65 Diamond Coated Forceps - Angled/Diamond Coated/14 cm



Diamond coating provides increased corrosion resistance and greater protection from wear and tear, thus extending the operational life of the forceps, as well as providing a stronger grip at the tips. All featured products are limited in quantity and may not be available in the future or part of the regular instrument line. While supplies last.

General Info

Product Number	11991-14
Category	Forceps & Tweezers > Smooth MedTech Line
EAN	4251909510220



Find 11991-14 on finescience.com

Technical Info

Shipping Weight	20 g	Net Weight	12 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Diamond coated Surgical Stainless Steel (Inox)		
Length	14 cm		
Tipdimensions	0.15 x 0.1 mm		
Tips	Standard		
Tipshape	Angled 30°		