



Halsey Needle Holder - Straight/Smooth/13cm/with Lock



General Info

Product Number	12001-13
Category	Needles, Threads & Holders > Needle Holders
EAN	4061706044442



Find 12001-13 on
finescience.com

Technical Info

Shipping Weight	30 g	Net Weight	26 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Clampinglength	16 mm		
Length	13 cm		
Feature	Smooth		
Lockable	Yes		
Scissors	No		
Serrated	No		
Tipwidth	1.9 mm		
Tipshape	Straight		