



Olsen-Hegar Needle Holder with Suture Cutter -
Straight/Serrated/12cm/with Lock



MOST
POPULAR



SERRATED

General Info

Product Number	12002-12
Category	Needles, Threads & Holders > Needle Holders Most Popular > Needles, Threats & Holders
EAN	4251909502768



Find 12002-12 on
finescience.com

Technical Info

Shipping Weight	20 g	Net Weight	20 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Clampinglength	8 mm		
Cuttingedge	10 mm		
Length	12 cm		
Lockable	Yes		
Scissors	Yes		
Serrated	Yes		
Tipwidth	1.5 mm		
Tipshape	Straight		