



Olsen-Hegar Needle Holder with Suture Cutter -
Straight/Serrated/14cm/with Lock



General Info

Product Number	12002-14
Category	Needles, Threads & Holders > Needle Holders Most Popular > Needles, Threats & Holders
EAN	4061706044930



Find 12002-14 on
finescience.com

Technical Info

Shipping Weight	30 g	Net Weight	27 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Clampinglength	8 mm		
Cuttingedge	15 mm		
Length	14 cm		
Lockable	Yes		
Scissors	Yes		
Serrated	Yes		
Tipwidth	2 mm		
Tipshape	Straight		