



Mayo-Hegar Needle Holder - Straight/Serrated/15cm/with Lock



Allgemein

|               |                                             |
|---------------|---------------------------------------------|
| Produktnummer | 12004-16                                    |
| Kategorie     | Needles, Threads & Holders > Needle Holders |
| EAN           | 4251909502782                               |



Finde 12004-16 auf [finescience.com](https://finescience.com)

Technische Informationen

|                     |                                 |              |      |
|---------------------|---------------------------------|--------------|------|
| Versandgewicht      | 45 g                            | Nettogewicht | 45 g |
| Paketgröße          | 1 Pc.                           | Maßeinheit   | EA   |
| Mindestbestellmenge | 1                               |              |      |
| Alloymaterial       | Surgical Stainless Steel (Inox) |              |      |
| Clampinglength      | 14 mm                           |              |      |
| Length              | 15 cm                           |              |      |
| Lockable            | Yes                             |              |      |
| Scissors            | No                              |              |      |
| Serrated            | Yes                             |              |      |
| Tipwidth            | 2 mm                            |              |      |
| Tipshape            | Straight                        |              |      |