



Slim Elongated Needle Holder - Tungsten Carbide/Straight/Serrated/15cm/with Lock



General Info

Product Number	12005-15
Category	Needles, Threads & Holders > Needle Holders
EAN	4251909502805



Find 12005-15 on finescience.com

Technical Info

Shipping Weight	30 g	Net Weight	24 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox) with tungsten carbide		
Clampinglength	9 mm		
Length	15 cm		
Feature	Tungsten Carbide jaws		
Lockable	Yes		
Scissors	No		
Serrated	Yes		
Tipwidth	1 mm		
Tipshape	Straight		