



Mathieu Needle Holder - Straight/Serrated/14cm/with Lock



General Info

Product Number	12010-14
Category	Needles, Threads & Holders > Needle Holders
EAN	4251909502812



Find 12010-14 on
[finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	40 g	Net Weight	35 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Clampinglength	14 mm		
Length	14 cm		
Feature	Extra narrow profile		
Lockable	Yes		
Scissors	No		
Serrated	Yes		
Tipwidth	2.3 mm		
Tipshape	Straight		