



Michel Suture Clip Applicator & Remover - For 11mm+ Clips



General Info

Product Number	12028-12
Category	Clamps & Clips > Michel System
EAN	4251909502874



Find 12028-12 on
finescience.com

Technical Info

Shipping Weight	35 g	Net Weight	31 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Length	13 cm		
Sterile	No		
Tipshape	Curved		