



Surgical Staple Kit - Handle/Remover/12 Non-Sterile Staple Cartridges



This reloadable skin stapler with a 360 swivel nose offers maximum visibility at the delivery point. The disposable staple cartridge includes a printed counter so you know how many staples you have left. Non-sterile, 316L stainless steel staples feature the patented Twin Towers design for quick and easy removal. The handle is reusable and autoclavable. Kit includes: 1 Staple Handle (No. 12030-01) 12 Non-Sterile Staple Cartridges (No. 12030-23) 1 Staple Remover (No. 12030-04) Kit components are also available separately.

General Info

Product Number	12030-20
Category	Clamps & Clips > FST System
EAN	4251909502911



Find 12030-20 on
finescience.com

Technical Info

Shipping Weight	555 g	Net Weight	530 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Autoclavesafe	Partial		
Sterile	No		