



Single Non-Sterile Staple Cartridge (1 Pack of 35 staples)



This disposable cartridge contains 35 staples and includes a printed counter so you know how many staples you have left. Non-sterile, 316L stainless steel staples for quick and easy handling. For use with Surgical Staple Handle (12030-01).

General Info

Product Number	12030-22
Category	Clamps & Clips > FST System
EAN	4251909502928



Find 12030-22 on
finescience.com

Technical Info

Shipping Weight	25 g	Net Weight	22 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	No		
Dimensionstext	6.5 x 4.8 mm		
Sterile	No		