



Reflex® 9 mm - Wound Clips (Pack of 100)



Suture Clips are easy and rapidly to apply with minimal trauma. The stainless steel clips are non-pyrogenic, have been associated with decreased wound infection rates and are less reactive than traditional suturing material. The applicator holds 20 clips which are preloaded on a rack. The clips will be supplied non-sterile. The applicator and remover are effortlessly to use and when the wound is healed, clips are easy to remove with the associated remover.

General Info

Product Number	12032-09
Category	Clamps & Clips > Reflex® System
EAN	4251909509101



Find 12032-09 on
finescience.com

Technical Info

Shipping Weight	25 g	Net Weight	23 g
Package Size	100 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Width	9 mm		
Sterile	No		