



Michel Suture Clips - 7.5 x 1.75mm

Made from 18/8 non-magnetic, acid resistant stainless steel, these clips are easily applied with a forceps-style applicator.



General Info

Product Number	12040-01
Category	Clamps & Clips > Michel System
EAN	4251909502942



Find 12040-01 on
finescience.com

Technical Info

Shipping Weight	95 g	Net Weight	91 g
Package Size	1000 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Dimensionstext	7.5 x 1.75 mm		
Sterile	No		