



Michel Suture Clips - 12 x 2.5mm

Made from 18/8 non-magnetic, acid resistant stainless steel, these clips are easily applied with a forceps-style applicator.



General Info

Product Number	12040-03
Category	Clamps & Clips > Michel System
EAN	4251909502966



Find 12040-03 on finescience.com

Technical Info

Shipping Weight	165 g	Net Weight	161 g
Package Size	1000 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Dimensionstext	12 x 2.5 mm		
Sterile	No		