



Kalt Suture Needles - Size 2



These Swiss made suture needles are uniform spring tempered, round, polished and fitted with a standard eye.

General Info

Product Number	12050-02
Category	Needles, Threads & Holders > Needles Most Popular > Needles, Threats & Holders
EAN	4251909502980



Find 12050-02 on
finescience.com

Technical Info

Shipping Weight	1 g	Net Weight	1 g
Package Size	12 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Stainless Steel (Inox)		
Autoclavesafe	Yes		
Needlechord	16 mm		
Needlecircle	3/8		
Needlediameter	0.35 mm		
Size	2		
Sterile	No		