




Microsurgical Needles with Suture Thread Attached - 8/0/0.045 mm



Allgemein

Produktnummer	12052-08
Kategorie	Needles, Threads & Holders > Needles with Thread
EAN	4251909511302



Finde 12052-08 auf
finescience.com

Technische Informationen

Versandgewicht	100 g	Nettogewicht	95 g
Paketgröße	12 Pc.	Maßeinheit	PA
Mindestbestellmenge	1		
Alloymaterial	Nylon		
Needlealloy	Stainless Steel		
Needlechord	29 mm		
Needlecircle	3/8		
Needlediameter	0.14 mm / 140 µm		
Size	8/0		
Sterile	Yes		
Suturediameter	0.045 mm		
Suturelength	15 cm		
Suturematerial	Nylon		