




Microsurgical Needles with Suture Thread Attached - 8/0/0.045 mm



General Info

Product Number	12052-08
Category	Needles, Threads & Holders > Needles with Thread
EAN	4251909511302



Find 12052-08 on
finescience.com

Technical Info

Shipping Weight	100 g	Net Weight	95 g
Package Size	12 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Nylon		
Autoclavesafe	Yes		
Needlealloy	Stainless Steel		
Needlechord	29 mm		
Needlecircle	3/8		
Needlediameter	0.14 mm / 140 µm		
Size	8/0		
Sterile	Yes		
Suturediameter	0.045 mm		
Suturelength	15 cm		
Suturematerial	Nylon		