




Microsurgical Needles with Suture Thread Attached - 9/0/0.035 mm



Allgemein

Produktnummer	12052-09
Kategorie	Needles, Threads & Holders > Needles with Thread
EAN	4251909511319



Technische Informationen

Versandgewicht	100 g	Nettogewicht	95 g
Paketgröße	12 Pc.	Maßeinheit	PA
Mindestbestellmenge	1		
Alloymaterial	Nylon		
Needlealloy	Stainless Steel		
Needlechord	24 mm		
Needlecircle	3/8		
Needlediameter	0.1 mm / 100 µm		
Size	9/0		
Sterile	Yes		
Suturediameter	0.035 mm		
Suturelength	15 cm		
Suturematerial	Nylon		