




Microsurgical Needles with Suture Thread Attached - 10/0/0.025 mm



Allgemein

Produktnummer	12052-10
Kategorie	Needles, Threads & Holders > Needles with Thread
EAN	4251909511326



Finde 12052-10 auf
finescience.com

Technische Informationen

Versandgewicht	100 g	Nettogewicht	95 g
Paketgröße	12 Pc.	Maßeinheit	PA
Mindestbestellmenge	1		
Alloymaterial	Nylon		
Autoclavesafe	Yes		
Needlealloy	Stainless Steel		
Needlechord	18 mm		
Needlecircle	3/8		
Needlediameter	0.07 mm / 70 µm		
Size	10/0		
Sterile	Yes		
Suturediameter	0.025 mm		
Suturelength	10 cm		
Suturematerial	Nylon		