




Microsurgical Needles with Suture Thread Attached - 11/0/0.015 mm



Allgemein

Produktnummer	12052-11
Kategorie	Needles, Threads & Holders > Needles with Thread
EAN	4251909511333



Finde 12052-11 auf finescience.com

Technische Informationen

Versandgewicht	100 g	Nettogewicht	95 g
Paketgröße	12 Pc.	Maßeinheit	PA
Mindestbestellmenge	1		
Alloymaterial	Nylon		
Needlealloy	Stainless Steel		
Needlechord	16 mm		
Needlecircle	3/8		
Needlediameter	0.05 mm / 50 µm		
Size	11/0		
Sterile	Yes		
Suturediameter	0.015 mm		
Suturelength	10 cm		
Suturematerial	Nylon		