




Microsurgical Needles with Suture Thread Attached - 11/0/0.015 mm



General Info

Product Number	12052-11
Category	Needles, Threads & Holders > Needles with Thread
EAN	4251909511333



Find 12052-11 on
finescience.com

Technical Info

Shipping Weight	100 g	Net Weight	95 g
Package Size	12 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Nylon		
Needlealloy	Stainless Steel		
Needlechord	16 mm		
Needlecircle	3/8		
Needlediameter	0.05 mm / 50 µm		
Size	11/0		
Sterile	Yes		
Suturediameter	0.015 mm		
Suturelength	10 cm		
Suturematerial	Nylon		