



Castroviejo Micro Needle Holder - Straight/Smooth/9cm/with Lock



Micro-pattern with slightly blunted tips. Cautionary note: Never use with needles that are too thick as this can easily damage the fine jaws. Use only micro needles of less than 0.2mm diameter.

General Info

Product Number	12060-01
Category	Needles, Threads & Holders > Needle Holders MedTech Line Most Popular > Needles, Threats & Holders
EAN	4251909503062



Find 12060-01 on
finescience.com

Technical Info

Shipping Weight	30 g	Net Weight	7 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Clampinglength	5 mm		
Length	9 cm		
Lockable	Yes		
Scissors	No		
Serrated	No		
Tipwidth	0.6 mm		
Tipshape	Straight		