



Round Handled Micro Needle Holder - Curved/Smooth/12.5cm/without Lock



Allgemein

Produktnummer	12074-12
Kategorie	Needles, Threads & Holders > Needle Holders
EAN	4251909511593



Finde 12074-12 auf finescience.com

Technische Informationen

Versandgewicht	47 g	Nettogewicht	13 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Clampinglength	5 mm		
Length	12.5 cm		
Lockable	No		
Scissors	No		
Serrated	No		
Tipwidth	0.6 mm		
Tipshape	Curved		