



Round Handled Micro Needle Holder with Suture Cutter - Curved/Smooth/14cm/without Lock



Cautionary note: Never use with needles that are too thick as this can easily damage the fine jaws. Use only micro needles of less than 0.2mm diameter.

General Info

Product Number	12075-14
Category	Needles, Threads & Holders > Needle Holders
EAN	4251909503116

Find 12075-14 on finescience.com

Technical Info

Shipping Weight	55 g	Net Weight	19 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Clampinglength	5 mm		
Cuttingedge	3.5 mm		
Length	14 cm		
Lockable	No		
Scissors	Yes		
Serrated	No		
Tipwidth	0.8 mm		
Tipshape	Curved		