



Round Handled Micro Needle Holder - Straight/Smooth/12.5cm/without Lock



Round handles allow easy fingertip adjustments. Cautionary note: Never use with needles that are too thick as this can easily damage the fine jaws. Use only micro needles of less than 0.2mm diameter.

General Info

Product Number	12076-12
Category	Needles, Threads & Holders > Needle Holders
EAN	4251909503123



Find 12076-12 on finescience.com

Technical Info

Shipping Weight	30 g	Net Weight	15 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Clampinglength	7 mm		
Length	12.5 cm		
Lockable	No		
Scissors	No		
Serrated	No		
Tipwidth	0.6 mm		
Tipshape	Straight		