



Kalt Needle Holder - Curved/Serrated/14cm/with Lock



The Kalt needle holder is characterized by its extraordinary design and good handling. The closure system closes softly and securely and can be easily opened again with your thumb. The needle is held securely when sewing thanks to various groove systems in the gripping surfaces (serration and cross-serration with blood channel). The needle holder is made of stainless steel and can be autoclaved. A straight and a curved variant is available.

General Info

Product Number	12081-14
Category	Needles, Threads & Holders > Needle Holders
EAN	4251909511753



Find 12081-14 on finescience.com

Technical Info

Shipping Weight	65 g	Net Weight	32 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Clampinglength	10 mm		
Length	14 cm		
Lockable	Yes		
Scissors	No		
Serrated	Yes		
Tipwidth	1 mm		
Tipshape	Curved		