



Halsey Needle Holder - Tungsten Carbide/Straight/Serrated/12cm/with Lock



MOST  
POPULAR



TUNGSTEN  
CARBIDE



SERRATED

General Info

|                |                                                                                          |
|----------------|------------------------------------------------------------------------------------------|
| Product Number | 12500-12                                                                                 |
| Category       | Needles, Threads & Holders > Needle Holders<br>Most Popular > Needles, Threats & Holders |
| EAN            | 4251909503147                                                                            |



Find 12500-12 on  
[finescience.com](https://finescience.com)

Technical Info

|                        |                                                       |                 |      |
|------------------------|-------------------------------------------------------|-----------------|------|
| Shipping Weight        | 20 g                                                  | Net Weight      | 17 g |
| Package Size           | 1 Pc.                                                 | Unit of Measure | EA   |
| Minimum Order Quantity | 1                                                     |                 |      |
| Alloymaterial          | Surgical Stainless Steel (Inox) with tungsten carbide |                 |      |
| Clampinglength         | 12.5 mm                                               |                 |      |
| Length                 | 12 cm                                                 |                 |      |
| Feature                | Tungsten Carbide jaws                                 |                 |      |
| Lockable               | Yes                                                   |                 |      |
| Scissors               | No                                                    |                 |      |
| Serrated               | Yes                                                   |                 |      |
| Tipwidth               | 1.4 mm                                                |                 |      |
| Tips                   | Cross Serrated Tungsten Carbide                       |                 |      |
| Tipshape               | Straight                                              |                 |      |