



Halsey Needle Holder - Tungsten Carbide/Straight/Serrated/13cm/with Lock



General Info

Product Number	12501-13
Category	Needles, Threads & Holders > Needle Holders
EAN	4251909503154



Find 12501-13 on finescience.com

Technical Info

Shipping Weight	30 g	Net Weight	28 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox) with tungsten carbide		
Clampinglength	15 mm		
Length	13 cm		
Feature	Tungsten Carbide jaws		
Lockable	Yes		
Scissors	No		
Serrated	Yes		
Tipwidth	2 mm		
Tips	Cross Serrated Tungsten Carbide		
Tipshape	Straight		