



Olsen-Hegar Needle Holder with Suture Cutter - Tungsten Carbide/Straight/Serrated/14cm/with Lock



General Info

Product Number	12502-14
Category	Needles, Threads & Holders > Needle Holders
EAN	4061706043957



Find 12502-14 on finescience.com

Technical Info

Shipping Weight	30 g	Net Weight	29 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox) with tungsten carbide		
Clampinglength	8 mm		
Cuttingedge	15 mm		
Length	14 cm		
Feature	Tungsten Carbide jaws		
Lockable	Yes		
Scissors	Yes		
Serrated	Yes		
Tipwidth	2 mm		
Tips	Cross Serrated Tungsten Carbide		
Tipshape	Straight		