



Crile-Wood Needle Holder - Tungsten Carbide/Straight/Serrated/15cm/with Lock



General Info

Product Number	12503-15
Category	Needles, Threads & Holders > Needle Holders
EAN	4251909503178



Find 12503-15 on [finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	30 g	Net Weight	29 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox) with tungsten carbide		
Clampinglength	14 mm		
Length	15 cm		
Lockable	Yes		
Scissors	No		
Serrated	Yes		
Tipwidth	2 mm		
Tips	Cross Serrated Tungsten Carbide		
Tipshape	Straight		