



Crile-Wood Needle Holder (Left-handed) - Tungsten Carbide/Straight/Serrated/15cm/with Lock/Left-Handed



Allgemein

Produktnummer	12504-15
Kategorie	Needles, Threads & Holders > Needle Holders
EAN	4251909503185



Finde 12504-15 auf finescience.com

Technische Informationen

Versandgewicht	35 g	Nettogewicht	29 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox) with tungsten carbide		
Clampinglength	14 mm		
Length	15 cm		
Islefthanded	Yes		
Lockable	Yes		
Scissors	No		
Serrated	Yes		
Tipwidth	2 mm		
Tips	Cross Serrated Tungsten Carbide		
Tipshape	Straight		