



Ryder Mini-Select Needle Holder - Tungsten Carbide Dust
Tips/Straight/15cm/with Lock



Tungsten Carbide Dust Tips give the clamping platform an extra secure grip.

This product is discontinued and has limited availability.

General Info

Product Number	12505-15
Category	Needles, Threads & Holders > Needle Holders
EAN	4251909503192



Find 12505-15 on
finescience.com

Technical Info

Shipping Weight	30 g	Net Weight	27 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox) with tungsten carbide		
Clampinglength	15 mm		
Length	15 cm		
Feature	Tungsten Carbide jaws		
Lockable	Yes		
Scissors	Yes		
Serrated	Yes		
Tipwidth	1.4 mm		
Tips	Cross Serrated Tungsten Carbide		
Tipshape	Straight		