



Mathieu Needle Holder - Tungsten Carbide/Straight/Serrated/14cm/with Lock



General Info

Product Number	12510-14
Category	Needles, Threads & Holders > Needle Holders
EAN	4251909503208



Find 12510-14 on finescience.com

Technical Info

Shipping Weight	35 g	Net Weight	33 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox) with tungsten carbide		
Clampinglength	14 mm		
Length	14 cm		
Feature	Tungsten Carbide jaws		
Lockable	Yes		
Scissors	No		
Serrated	Yes		
Tipwidth	2.3 mm		
Tips	Cross Serrated Tungsten Carbide		
Tipshape	Straight		