



Olsen-Hegar Needle Holder with Suture Cutter - Tungsten Carbide/Straight/Serrated/12cm/with Lock/Left-Handed



Allgemein

Produktnummer	12512-12
Kategorie	Needles, Threads & Holders > Needle Holders
EAN	4251909503215



Finde 12512-12 auf finescience.com

Technische Informationen

Versandgewicht	20 g	Nettogewicht	16 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox) with tungsten carbide		
Clampinglength	8 mm		
Cuttingedge	10 mm		
Length	12 cm		
Feature	Tungsten Carbide jaws		
Islefthanded	Yes		
Lockable	Yes		
Scissors	Yes		
Serrated	Yes		
Tipwidth	1.5 mm		
Tips	Cross Serrated Tungsten Carbide		
Tipshape	Straight		