



Castroviejo Needle Holder - Tungsten Carbide/Straight/Serrated/14cm/with Lock



General Info

Product Number	12565-14
Category	Needles, Threads & Holders > Needle Holders MedTech Line
EAN	4251909503222



Find 12565-14 on
finescience.com

Technical Info

Shipping Weight	55 g	Net Weight	20 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox) with tungsten carbide		
Clampinglength	10 mm		
Length	14 cm		
Feature	Tungsten Carbide		
Lockable	Yes		
Scissors	No		
Serrated	Yes		
Tipwidth	1.5 mm		
Tips	Cross Serrated Tungsten Carbide		
Tipshape	Straight		