



Evo Needle Holder - Straight/Smooth/15cm/without Lock



An evolutionary series of ergonomic instruments that are lighter than titanium but still have all the benefits of stainless steel. The fully fenestrated spring body is designed with a stability rod to provide balance in your hand. Cautionary note: Never use with needles that are too thick as this can easily damage the fine jaws. Use only micro needles of less than 0.2mm diameter.

This product is discontinued and has limited availability.

General Info

Product Number	12800-00
Category	Needles, Threads & Holders > Needle Holders
EAN	4251909509125



Find 12800-00 on finescience.com

Technical Info

Shipping Weight	50 g	Net Weight	12 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Clampinglength	5 mm		
Length	15 cm		
Lockable	No		
Scissors	No		
Serrated	No		
Tipwidth	0.6 mm		
Tipshape	Straight		