



Round-Handled Micro/Macro Needle Holder -
Curved/Smooth/13.5cm/without Lock



Our versatile Round Handled Needle Holders of the Exclusive Line are designed for both microscopic and macroscopic tasks. Finely tapered tips securely grip small needles up to .2mm in diameter and also ensure the same level of precision with larger needles using the sturdy rear section. The rounded handle enables swift and precise movement with unmatched control.

General Info

Product Number	12875-13
Category	Exclusive Line Needles, Threads & Holders > Needle Holders MedTech Line
EAN	4251909511975



Find 12875-13 on
finescience.com

Technical Info

Shipping Weight	55 g	Net Weight	19 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Clampinglength	9 mm		
Length	13.5 cm		
Lockable	No		
Scissors	No		
Serrated	No		
Tipwidth	0.4 mm		
Tipshape	Curved		