



Hartmann Hemostat - Serrated/Straight/10cm



Allgemein

Produktnummer	13002-10
Kategorie	Clamps & Clips > Hemostats Most Popular > Clamps & Clips
EAN	4251909503239



Finde 13002-10 auf finescience.com

Technische Informationen

Versandgewicht	15 g	Nettogewicht	13 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Clampinglength	18 mm		
Length	10 cm		
Tipwidth	1 mm		
Tips	Serrated		
Tipshape	Straight		